


DopUS-Net: Quality-Aware Robotic Ultrasound Imaging Based on Doppler Signal

Zhongliang Jiang*, Felix Duellmer*, and Nassir Navab, *Fellow, IEEE*

Affiliation: Computer Aided Medical Procedures, Technical University of Munich, Munich, Germany

*Contributed equally



ICRA2024
YOKOHAMA | JAPAN



Video



Code



Paper

Motivation

Peripheral Artery Disease (PAD) is associated with:

- Increased risk of cardiovascular morbidity and mortality
- Reduced functional capacity


It is estimated that >200 million people have PAD worldwide [1]

Regular screening procedures are recommended, but a lack of clinical experts and expensive equipment limit their application.

→ We propose a robust quality-aware robotic Ultrasound screening using Doppler information for autonomous peripheral artery segmentation and 3D reconstruction.

DopUS-Net

Segmentation



Two Encoders: enable the optimal use of the Doppler effect:

- Top encoder uses a 2-channel input of B-Mode and Doppler grey-scale image for enhanced structural accuracy
- Bottom encoder uses only Doppler images to act as a region proposal module

Temporal Continuity: exploits the anatomical continuity of vessels

- A convolutional gated recurrent unit (ConvGRU) uses a hidden state to keep information from previous frames in memory

Quantitative Evaluation

TABLE II
COMPARISON OF THE RESULTS

| Network | Top Encoder | Bottom Encoder | # Parameters | 0 | 1 | 2 | Patient 3 | 4 | 5 | 6 | Dice Score Mean (SD) |
|--------------------------|-------------|----------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|----------------------|
| U-Net | B | - | 0.6 M | 0.61 | 0.51 | 0.45 | 0.40 | 0.45 | 0.29 | 0.38 | 0.44 (0.09) |
| U-Net | BD | - | 0.6 M | 0.83 | 0.59 | 0.54 | 0.52 | 0.49 | 0.39 | 0.45 | 0.55 (0.13) |
| DopUS-Net ⁽⁰⁾ | BD | D | 1.3 M | 0.83 | 0.61 | 0.58 | 0.60 | 0.54 | 0.43 | 0.47 | 0.58 (0.12) |
| VesNet | BD-RNN | - | 2.6 M | 0.80 | 0.62 | 0.60 | 0.43 | 0.33 | 0.42 | 0.43 | 0.52 (0.15) |
| VesNet+ | BD-RNN | - | 6.3 M | 0.80 | 0.62 | 0.68 | 0.58 | 0.58 | 0.47 | 0.53 | 0.61 (0.10) |
| U-Net | BD-RNN | - | 3.0 M | 0.86 | 0.61 | 0.71 | 0.66 | 0.68 | 0.60 | 0.50 | 0.66 (0.10) |
| U-Net+ | BD-RNN | - | 6.5 M | 0.78 | 0.58 | 0.37 | 0.40 | 0.36 | 0.32 | 0.41 | 0.46 (0.15) |
| DopUS-Net ⁽¹⁾ | B-RNN | D-RNN* | 6.2 M | 0.78 | 0.54 | 0.45 | 0.56 | 0.55 | 0.27 | 0.18 | 0.48 (0.19) |
| DopUS-Net ⁽²⁾ | BD | D-RNN | 3.7 M | 0.87 | 0.69 | 0.76 | 0.76 | 0.72 | 0.56 | 0.61 | 0.71 (0.10) |
| DopUS-Net ⁽³⁾ | B-RNN | D-RNN | 6.1 M | 0.88 | 0.71 | 0.79 | 0.76 | 0.75 | 0.57 | 0.61 | 0.72 (0.10) |
| DopUS-Net ⁽⁴⁾ | BD-RNN | D-RNN | 6.1 M | 0.88 | 0.69 | 0.79 | 0.78 | 0.76 | 0.62 | 0.60 | 0.73 (0.09) |

*Nomenclature: B: B-Mode, D: Doppler, RNN: convGRU module, +: increased parameters, DopUS-Net⁽⁰⁾: specific DopUS-Net version, *: additional skip connections from the bottom encoder to the decoder.

Performance: Compared with U-Net [2] and VesNet [3], DopUS-Net outperforms the other networks in terms of Dice score.

Ablations & Observations:


- Leveraging an RNN module enhances the performance. Pulsation of blood flow leads to unstable Doppler signal, which can be alleviated through temporal information.
- Using a separate encoder for the Doppler signal enhances the performance

Conclusion

We introduced a quality-aware robotic Ultrasound screening using Doppler information for autonomous peripheral artery segmentation and 3D reconstruction:

- DopUS-Net: novel segmentation network leveraging Doppler and continuity information for superior segmentation performance
- Quality-aware scanning: Doppler re-identification procedure for robust reconstruction performance using a closed-loop control scheme



System Setup



Doppler Re-Identification

Re-Identification

Doppler effect in the forearm is most prominent in the vessels.



Quality-Awareness: This module constantly checks if the predicted segmentation mask aligns with the Doppler signal. In case of an insufficient overlap a re-identification procedure is triggered.


Re-identification procedure:

Goal: Re-align segmentation mask with Doppler signal

Procedure:

- Rotate in 5° steps in the out-of-plane direction: [-10°, -5°, 0°, 5°, 10°]
- After every fan motion (rotation) check for re-alignment and artery re-identification
- If no step leads to a re-identification continue scan

Qualitative Evaluation



Quantitative Evaluation

Performance: Re-identification procedure enhances overall Dice score and IoU due to the quality-awareness module.

TABLE III
SEGMENTATION RESULTS OVER VOLUME [MEAN (SD)]

| Re-Identification | Dice Score | IoU |
|-------------------|-------------|-------------|
| Enabled | 0.86 (0.14) | 0.78 (0.16) |
| Disabled | 0.54 (0.39) | 0.47 (0.37) |

References

- [1] Shu, J., & Santulli, G. (2018). Update on peripheral artery disease: Epidemiology and evidence-based facts. *Atherosclerosis*, 275, 379-381.
- [2] Ronneberger, O., Fischer, P., & Brox, T. (2015). U-net: Convolutional networks for biomedical image segmentation. In *Medical image computing and computer-assisted intervention—MICCAI 2015: 18th international conference, Munich, Germany, October 5-9, 2015, proceedings, part III* 18 (pp. 234-241). Springer International Publishing.
- [3] Jiang, B., Chen, A., Bharat, S., & Zheng, M. (2021). Automatic ultrasound vessel segmentation with deep spatiotemporal context learning. In *Simplifying Medical Ultrasound: Second International Workshop, ASMUS 2021, Held in Conjunction with MICCAI 2021, Strasbourg, France, September 27, 2021, Proceedings 2* (pp. 3-13). Springer International Publishing.